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OPERATING INSTRUCTIONS

3.1 Controls, Terminals and Indicators

a. FRONT PANEL

GATE ON or PRINT Indicator : Lights to indicate an open gate or interrogation in process.

OVERRANGE Indicators : Lights to indicate the scaler has been overranged.

MASTER/SLAVE Switch : In MASTER position the scaler is controlled either by the front-panel controls or by signals applied to the rear-panel connectors. A gate output is provided to control the slave modules.

In SLAVE position the scaler is controlled by front-panel controls and by external signals applied to the rear-panel connectors by the master of the multiscaler system.

NORM/INTEG Switch : In NORM position reset pulses are enabled. In INTEG position the counting is cumulative.

PRESET Switches : Two rotary switches for selection of the preset number. In the OFF position the counting can be stopped only by the STOP push-button, or externally.

RECY/NORM Switch : Enables recycle mode.

PRESET TIME/COUNT Switch : Selects either preset time or count time mode.

DISPLAY TIME/COUNT Switch : Selects either the time or count for display.

START PUSHBUTTON : Starts the counting.

NORM/DISP.TEST/TEST Switch : In the DISP.TEST position the display indicates 88888 effecting the count. In the TEST position a 1 MHz signal is applied to the counter scaler input and is counted for the duration of the gate-pulse.

STOP PUSHBUTTON : Stops the counting.

SENSE. CONTROLS : Adjusts the threshold of the input pulse from +12 volts to -12 volts.

DISPLAY TIME Control : In RECY. mode, adjusts the time between cycles, from 0.5 to 10 seconds. However, it will not restart until the printing is complete.

INPUT BNC : For connection of input signal.

b. REAR PANEL

PRINTER OUT Connection : : 50-pin connector for output to printer or via SPC-1 to TTY. Contains seven BCD digits, print command and other control signals used in a multiscaler system.

ON/OFF Switch : In the OFF position the scaler is excluded from a multiscaler system by turning off the 5-volt power.

GATE IN BNC : Control starts and stops counting.

GATE OUT BNC : Output pulse during counting.

START IN BNC : Starts counting on minimum of 1 μ sec pulse.

STOP IN BNC : Stops counting on minimum of 1 μ sec pulse.

OVERFLOW BNC : Pulse output when the capacity of the scaler has been exceeded.

CLOCK INT/EXT Switch : In EXT position clocks are supplied to CLOCK IN/OUT BNC. In INT position, internal clocks are used and are made available via the CLOCK IN/OUT BNC.

CLOCK IN/OUT BNC : Depending on position of CLOCK INT/EXT switch, accepts external clocks or supplies internal clocks.

POWER CONNECTOR : Standard connector for NIM power supply such as ELSCINT BPS-N-3.

System Connector:

For interconnection to a multiscaler system or to Terminal Box TBS-1.

3.2 Operation

3.2.1 Connect the STP-N-75A to a standard NIM power supply such as ELSCINT BPS-N-3.

3.2.2 Input circuit adjustment:

- a) Connect input signal to INPUT BNC on Front Panel, and place ON/OFF switch on rear-panel to the ON position.
- b) Place DISPLAY COUNT/TIME Switch to the COUNT position.
- c) Depress START pushbutton, GATE ON or PRINT indicator should light.
- d) Slowly rotate SENS. Control clockwise until counting is indicated on the display.
- e) Stop counting by depressing STOP pushbutton. The input circuit of the STP-N-75A is now ready for counting.

3.2.3 Master Mode operation:

- a) Place MASTER/SLAVE switch to the MASTER position.
- b) Connect SLAVE modules , if any, to the rear-panel connector .
- c) For cumulative counting of the input signal place NORM/INTEG switch to the INTEG position.
- d) Connect external printer, if used, to the PRINTER OUT connector on the rear panel. If used with TTY page printer, connect Parallel to Serial Control SPC-1, to the rear cable.
- e) If the preset capability is to be used, set the desired preset number, for either time or count depending on position of PRESET TIME/COUNT switch, with the two controls on the front panel.
- f) Depress START pushbutton to start the counting. If the preset capability is used, counting will stop automatically when the preset number is reached. Otherwise, stop counting by depressing the STOP pushbutton.